

## Product datasheet for **TS430923P5**

### CENPS-CORT CytoSection

#### Product data:

|                                              |                                                                                                                                                                                                                                                                                                                                   |
|----------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b>                         | CytoSections                                                                                                                                                                                                                                                                                                                      |
| <b>Description:</b>                          | Transient overexpression of CENPS-CORT, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack                                                                                                                                                                                    |
| <b>Species:</b>                              | Human                                                                                                                                                                                                                                                                                                                             |
| <b>Expression Host:</b>                      | HEK293T                                                                                                                                                                                                                                                                                                                           |
| <b>Expression cDNA Clone or AA Sequence:</b> | TrueORF Clone RC230923                                                                                                                                                                                                                                                                                                            |
| <b>Tag:</b>                                  | C-MYC/DDK                                                                                                                                                                                                                                                                                                                         |
| <b>Detection Antibodies:</b>                 | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)                                                                                                                                                                                                                                                |
| <b>ACCN:</b>                                 | <u>NM_199006</u> , <u>NP_950171</u>                                                                                                                                                                                                                                                                                               |
| <b>Synonyms:</b>                             | APITD1; APITD1-CORT; CENP-S; CENPS; FAAP16; MHF1                                                                                                                                                                                                                                                                                  |
| <b>Storage:</b>                              | RT                                                                                                                                                                                                                                                                                                                                |
| <b>Stability:</b>                            | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.                                                                                                                                                                                                                           |
| <b>Preparation:</b>                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| <b>Note:</b>                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.                                                                                                                                                                                                                             |
| <b>RefSeq:</b>                               | <u>NP_950171</u>                                                                                                                                                                                                                                                                                                                  |
| <b>Locus ID:</b>                             | 100526739                                                                                                                                                                                                                                                                                                                         |
| <b>Cytogenetics:</b>                         | 1p36.22                                                                                                                                                                                                                                                                                                                           |

