

Product datasheet for **TS431279P5**

Sterol carrier protein 2 (SCP2) CytoSection

Product data:

| | |
|---------------------------------------|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of SCP2, transcript variant 6, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC231279 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| Target Detection Antibodies: | Sterol carrier protein 2 (SCP2) Mouse Monoclonal Antibody [Clone ID: OT11D2] (TA507248) |
| ACCN: | NM_001193600 , NP_001180529 |
| Synonyms: | NLTP; NSL-TP; SCOX; SCP-2; SCP-CHI; SCP-X; SCPX |
| Storage: | Room Temperature |
| Stability: | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | NP_001180529 |
| Locus ID: | 6342 |
| Cytogenetics: | 1p32.3 |
| Protein Pathways: | Metabolic pathways, PPAR signaling pathway, Primary bile acid biosynthesis |



[View online »](#)