

Product datasheet for **TS431376P5**

Dystrobrevin alpha (DTNA) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DTNA, transcript variant 11, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC231376
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Dystrobrevin alpha (DTNA) Mouse Monoclonal Antibody [Clone ID: OT11B2] (TA502184)
ACCN:	NM_001198939 , NP_001185868
Synonyms:	D18S892E; DRP3; DTN; DTN-A; LVNC1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001185868
Locus ID:	1837
Cytogenetics:	18q12.1
Protein Families:	Druggable Genome


[View online »](#)