

## Product datasheet for **TS431920**

### p21 (CDKN1A) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CDKN1A, transcript variant 5, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC231920
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	p21 (CDKN1A) Mouse Monoclonal Antibody [Clone ID: OTI5G1] (TA808128)
ACCN:	<u><a href="#">NM_001220777</a></u> , <u><a href="#">NP_001207706</a></u>
Synonyms:	CAP20; CDKN1; CIP1; MDA-6; P21; p21CIP1; SDI1; WAF1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001207706</a></u>
Locus ID:	1026
Cytogenetics:	6p21.2
Protein Families:	Druggable Genome
Protein Pathways:	Bladder cancer, Cell cycle, Chronic myeloid leukemia, ErbB signaling pathway, Glioma, Melanoma, p53 signaling pathway, Pathways in cancer, Prostate cancer



[View online »](#)