

## Product datasheet for **TS432013P5**

### Apc10 (ANAPC10) CytoSection

#### Product data:

|                                       |                                                                                                                                                                                                                                                                                                                                   |
|---------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:                         | CytoSections                                                                                                                                                                                                                                                                                                                      |
| Description:                          | Transient overexpression of ANAPC10, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack                                                                                                                                                                                       |
| Species:                              | Human                                                                                                                                                                                                                                                                                                                             |
| Expression Host:                      | HEK293T                                                                                                                                                                                                                                                                                                                           |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC232013                                                                                                                                                                                                                                                                                                            |
| Tag:                                  | C-MYC/DDK                                                                                                                                                                                                                                                                                                                         |
| Detection Antibodies:                 | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)                                                                                                                                                                                                                                                |
| ACCN:                                 | <u>NM_001256707</u> , <u>NP_001243636</u>                                                                                                                                                                                                                                                                                         |
| Synonyms:                             | APC10; DOC1                                                                                                                                                                                                                                                                                                                       |
| Storage:                              | Room Temperature                                                                                                                                                                                                                                                                                                                  |
| Stability:                            | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.                                                                                                                                                                                                                           |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.                                                                                                                                                                                                                             |
| RefSeq:                               | <u>NP_001243636</u>                                                                                                                                                                                                                                                                                                               |
| Locus ID:                             | 10393                                                                                                                                                                                                                                                                                                                             |
| Cytogenetics:                         | 4q31.21                                                                                                                                                                                                                                                                                                                           |
| Protein Families:                     | Druggable Genome                                                                                                                                                                                                                                                                                                                  |
| Protein Pathways:                     | Cell cycle, Oocyte meiosis, Progesterone-mediated oocyte maturation, Ubiquitin mediated proteolysis                                                                                                                                                                                                                               |



[View online »](#)