

Product datasheet for **TS432026P5**

CYBC1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CYBC1, transcript variant 7, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC232026
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CYBC1 Mouse Monoclonal Antibody [Clone ID: OTI3B5] (TA506979)
ACCN:	NM_001193657 , NP_001180586
Synonyms:	C17orf62; CGD5; Eros
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001180586
Locus ID:	79415
Cytogenetics:	17q25.3
Protein Families:	Transmembrane



[View online »](#)