

Product datasheet for **TS432315P5**

HLA-DQB1 CytoSection

Product data:

| | |
|---------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CytoSections |
| Description: | Transient overexpression of HLA-DQB1, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC232315 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| ACCN: | <u>NM_001243962</u> , <u>NP_001230891</u> |
| Synonyms: | CELIAC1; HLA-DQB; IDDM1 |
| Storage: | Room Temperature |
| Stability: | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | <u>NP_001230891</u> |
| Locus ID: | 3119 |
| Cytogenetics: | 6p21.32 |
| Protein Families: | Transmembrane |
| Protein Pathways: | Allograft rejection, Antigen processing and presentation, Asthma, Autoimmune thyroid disease, Cell adhesion molecules (CAMs), Graft-versus-host disease, Systemic lupus erythematosus, Type I diabetes mellitus, Viral myocarditis |



[View online »](#)