

Product datasheet for **TS432332P5**

ANKS1B CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ANKS1B, transcript variant 11, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC232332
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001204080</u> , <u>NP_001191009</u>
Synonyms:	AIDA; AIDA-1; ANKS2; cajalin-2; EB-1; EB1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001191009</u>
Locus ID:	56899
Cytogenetics:	12q23.1



[View online »](#)