

## Product datasheet for **TS432441**

### Asialoglycoprotein Receptor 2 (ASGR2) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ASGR2, transcript variant 4, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC232441
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Asialoglycoprotein Receptor 2 (ASGR2) Mouse Monoclonal Antibody [Clone ID: OTI2A12] (TA505099)
ACCN:	<u><a href="#">NM_001201352</a></u> , <u><a href="#">NP_001188281</a></u>
Synonyms:	ASGP-R2; ASGPR2; CLEC4H2; HBXBP; HL-2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001188281</a></u>
Locus ID:	433
Cytogenetics:	17p13.1
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)