

Product datasheet for **TS433222P5**

BCKDH kinase (BCKDK) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of BCKDK, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC233222
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	BCKDH kinase (BCKDK) Mouse Monoclonal Antibody [Clone ID: OT111C9] (TA500685)
ACCN:	NM_001271926 , NP_001258855
Synonyms:	BCKDKD; BDK
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001258855
Locus ID:	10295
Cytogenetics:	16p11.2
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)