

Product datasheet for **TS433980P5**

TRK fused gene (TFG) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TFG, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC233980
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	TRK fused gene (TFG) Mouse Monoclonal Antibody [Clone ID: OTI2C3] (TA503858)
ACCN:	NM_001195478 , NP_001182407
Synonyms:	HMSNP; SPG57; TF6; TRKT3
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001182407
Locus ID:	10342
Cytogenetics:	3q12.2
Protein Pathways:	Pathways in cancer, Thyroid cancer



[View online »](#)