

Product datasheet for **TS434032P5**

alpha Tubulin (TUBA1A) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TUBA1A, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC234032
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	alpha Tubulin (TUBA1A) Mouse Monoclonal Antibody [Clone ID: OTI2C8] (TA506605)
ACCN:	<u>NM_001270400</u> , <u>NP_001257329</u>
Synonyms:	B-ALPHA-1; LIS3; TUBA3
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001257329</u>
Locus ID:	7846
Cytogenetics:	12q13.12
Protein Families:	Druggable Genome
Protein Pathways:	Gap junction, Pathogenic Escherichia coli infection



[View online »](#)