

## Product datasheet for **TS434231**

### Cytokeratin 8 (KRT8) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of KRT8, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC234231
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Cytokeratin 8 (KRT8) Mouse Monoclonal Antibody [Clone ID: OT11B12] (TA500021)
ACCN:	<a href="#">NM_001256293</a> , <a href="#">NP_001243222</a>
Synonyms:	CARD2; CK-8; CK8; CYK8; K2C8; K8; KO
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001243222</a>
Locus ID:	3856
Cytogenetics:	12q13.13
Protein Families:	Druggable Genome



[View online »](#)