

Product datasheet for **TS434557P5**

Myoneurin (MYNN) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MYNN, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC234557
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Myoneurin (MYNN) Mouse Monoclonal Antibody [Clone ID: OT11F10] (TA800192)
ACCN:	NM_001185118 , NP_001172047
Synonyms:	OSZF; SBBIZ1; ZBTB31; ZNF902
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001172047
Locus ID:	55892
Cytogenetics:	3q26.2
Protein Families:	Transcription Factors



[View online »](#)