

Product datasheet for **TS434774P5**

Miz1 (ZBTB17) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ZBTB17, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC234774
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Miz1 (ZBTB17) Mouse Monoclonal Antibody [Clone ID: OTI7D8] (TA811343)
ACCN:	NM_001242884 , NP_001229813
Synonyms:	MIZ-1; pHZ-67; ZNF60; ZNF151
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001229813
Locus ID:	7709
Cytogenetics:	1p36.13
Protein Families:	Transcription Factors
Protein Pathways:	Cell cycle



[View online »](#)