

## Product datasheet for **TS435094P5**

### Myosin Phosphatase (PPP1R12A) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PPP1R12A, transcript variant 5, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC235094
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001244992</u> , <u>NP_001231921</u>
Synonyms:	GUBS; M130; MBS; MYPT1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001231921</u>
Locus ID:	4659
Cytogenetics:	12q21.2-q21.31
Protein Families:	Druggable Genome
Protein Pathways:	Focal adhesion, Long-term potentiation, Regulation of actin cytoskeleton, Vascular smooth muscle contraction



[View online »](#)