

Product datasheet for **TS435745P5**

Calcipressin 1 (RCAN1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RCAN1, transcript variant 7, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC235745
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Calcipressin 1 (RCAN1) Mouse Monoclonal Antibody [Clone ID: OTI10C6] (TA810470)
ACCN:	NM_001285393 , NP_001272322
Synonyms:	ADAPT78; CSP1; DSC1; DSCR1; MCIP1; RCN1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001272322
Locus ID:	1827
Cytogenetics:	21q22.12
Protein Families:	Transcription Factors



[View online »](#)