

## Product datasheet for **TS435839P5**

### **POLR1H CytoSection**

#### **Product data:**

<b>Product Type:</b>	CytoSections
<b>Description:</b>	Transient overexpression of POLR1H, transcript variant d, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	TrueORF Clone RC235839
<b>Tag:</b>	C-MYC/DDK
<b>Detection Antibodies:</b>	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
<b>Target Detection Antibodies:</b>	POLR1H Mouse Monoclonal Antibody [Clone ID: OT11G2] (TA505067)
<b>ACCN:</b>	<a href="#">NM_001278786</a> , <a href="#">NP_001265715</a>
<b>Synonyms:</b>	A12.2; HTEX-6; HTEX6; hZR14; Rpa12; tctex-6; TCTEX6; TEX6; ZNRD1; ZR14
<b>Storage:</b>	Room Temperature
<b>Stability:</b>	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
<b>Preparation:</b>	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
<b>Note:</b>	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
<b>RefSeq:</b>	<a href="#">NP_001265715</a>
<b>Locus ID:</b>	30834
<b>Cytogenetics:</b>	6p22.1
<b>Protein Families:</b>	Transcription Factors
<b>Protein Pathways:</b>	Metabolic pathways, Purine metabolism, Pyrimidine metabolism, RNA polymerase



[View online »](#)