

Product datasheet for **TS436269P5**

MASA (ENOPH1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ENOPH1, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC236269
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	MASA (ENOPH1) Mouse Monoclonal Antibody [Clone ID: OT16B6] (TA810832)
ACCN:	NM_001292017 , NP_001278946
Synonyms:	E1; MASA; MST145; mtnC
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001278946
Locus ID:	58478
Cytogenetics:	4q21.22
Protein Pathways:	Cysteine and methionine metabolism



[View online »](#)