

## Product datasheet for **TS437326P5**

### OR8G2P CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of OR8G2P in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC237326
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<a href="#">NM_001291438</a> , <a href="#">NP_001278367</a>
Synonyms:	HSTPCR120; OR8G2; OR8G4; ORL206; ORL486; TPCR120
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 $\mu$ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001278367</a>
Locus ID:	26492
Cytogenetics:	11q24.2
Protein Families:	Druggable Genome
Protein Pathways:	Olfactory transduction



[View online »](#)