

Product datasheet for **TS437407P5**

Cyclin Y (CCNY) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CCNY, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC237407
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Cyclin Y (CCNY) Mouse Monoclonal Antibody [Clone ID: OT112F10] (TA808174)
ACCN:	NM_001282852 , NP_001269781
Synonyms:	C10orf9; CBCP1; CCNX; CFP1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001269781
Locus ID:	219771
Cytogenetics:	10p11.21



[View online »](#)