

Product datasheet for **TS437408**

Apolipoprotein E (APOE) CytoSection

Product data:

| | |
|---------------------------------------|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of APOE, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC237408 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| Target Detection Antibodies: | Apolipoprotein E (APOE) Mouse Monoclonal Antibody [Clone ID: OT13B9] (TA805358) |
| ACCN: | NM_001302689 , NP_001289618 |
| Synonyms: | AD2; APO-E; ApoE4; LDLCQ5; LPG |
| Storage: | Room Temperature |
| Stability: | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | NP_001289618 |
| Locus ID: | 348 |
| Cytogenetics: | 19q13.32 |
| Protein Families: | Adult stem cells, Druggable Genome, Secreted Protein, Stem cell - Pluripotency |
| Protein Pathways: | Alzheimer's disease |



[View online »](#)