

## Product datasheet for **TS437489P5**

### GRAP2 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of GRAP2, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC237489
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	GRAP2 Mouse Monoclonal Antibody [Clone ID: OT11G2] (TA502923)
ACCN:	<a href="#">NM_001291825</a> , <a href="#">NP_001278754</a>
Synonyms:	GADS; GRAP-2; GRB2L; GRBLG; GrbX; Grf40; GRID; GRPL; Mona; P38
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001278754</a>
Locus ID:	9402
Cytogenetics:	22q13.1
Protein Families:	Druggable Genome
Protein Pathways:	T cell receptor signaling pathway



[View online »](#)