

## Product datasheet for **TS437567P5**

### G protein alpha Inhibitor 2 (GNAI2) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of GNAI2, transcript variant 6, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC237567
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_001282620</a></u> , <u><a href="#">NP_001269549</a></u>
Synonyms:	GIP; GNAI2B; H_LUCA15.1; H_LUCA16.1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001269549</a></u>
Locus ID:	2771
Cytogenetics:	3p21.31
Protein Families:	Druggable Genome
Protein Pathways:	Axon guidance, Chemokine signaling pathway, Gap junction, Leukocyte transendothelial migration, Long-term depression, Melanogenesis, Progesterone-mediated oocyte maturation, Tight junction



[View online »](#)