

## Product datasheet for **TS437651P5**

### ALDH3B1 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ALDH3B1, transcript variant 4, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC237651
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	ALDH3B1 Mouse Monoclonal Antibody [Clone ID: OTI2F6] (TA811463)
ACCN:	<a href="#">NM_001290058</a> , <a href="#">NP_001276987</a>
Synonyms:	ALDH4; ALDH7
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001276987</a>
Locus ID:	221
Cytogenetics:	11q13.2
Protein Families:	Druggable Genome



[View online »](#)

**Protein Pathways:**

Drug metabolism - cytochrome P450, Glycolysis / Gluconeogenesis, Histidine metabolism, Metabolic pathways, Metabolism of xenobiotics by cytochrome P450, Phenylalanine metabolism, Tyrosine metabolism