

Product datasheet for **TS437657P5**

NAT1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NAT1, transcript variant 10, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC237657
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	NAT1 Mouse Monoclonal Antibody [Clone ID: OTI6D7] (TA810719)
ACCN:	<u>NM_001291962</u> , <u>NP_001278891</u>
Synonyms:	AAC1; MNAT; NAT-1; NATI
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001278891</u>
Locus ID:	9
Cytogenetics:	8p22
Protein Pathways:	Caffeine metabolism, Drug metabolism - other enzymes, Metabolic pathways



[View online »](#)