

Product datasheet for **TS438346P5**

CDC2L6 (CDK19) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CDK19, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC238346
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001300960</u> , <u>NP_001287889</u>
Synonyms:	bA346C16.3; CDC2L6; CDK11; DEE87; EIEE87
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001287889</u>
Locus ID:	23097
Cytogenetics:	6q21
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)