

Product datasheet for **TS438567P5**

EAAT1 (SLC1A3) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SLC1A3, transcript variant GLAST1a, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC238567
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_001289939 , NP_001276868
Synonyms:	EA6; EAAT1; GLAST; GLAST1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001276868
Locus ID:	6507
Cytogenetics:	5p13.2
Protein Families:	Transmembrane



[View online »](#)