

## Product datasheet for **TS438803P5**

### Stromal interaction molecule 1 (STIM1) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of STIM1, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC238803
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_001277962</a></u> , <u><a href="#">NP_001264891</a></u>
Synonyms:	D11S4896E; GOK; IMD10; STRMK; TAM; TAM1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001264891</a></u>
Locus ID:	6786
Cytogenetics:	11p15.4
Protein Families:	Transmembrane



[View online »](#)