

## Product datasheet for **TS438922P5**

### MSF (SEPT9) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SEPTIN9, transcript variant 8, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC238922
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	MSF (SEPT9) Mouse Monoclonal Antibody [Clone ID: OT113B1] (TA812821)
ACCN:	<a href="#">NM_001293695</a> , <a href="#">NP_001280624</a>
Synonyms:	AF17q25; MSF; MSF1; NAPB; PNUTL4; SEPT9; SeptD1; SINT1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001280624</a>
Locus ID:	10801
Cytogenetics:	17q25.3
Protein Families:	Druggable Genome



[View online »](#)