

Product datasheet for **TS439049P5**

RPS6KC1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RPS6KC1, transcript variant 5, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC239049
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	RPS6KC1 Mouse Monoclonal Antibody [Clone ID: OTI10F11] (TA808825)
ACCN:	NM_001287220 , NP_001274149
Synonyms:	humS6PKh1; RPK118; RSKL1; S6K-delta-1; S6PKh1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001274149
Locus ID:	26750
Cytogenetics:	1q32.3
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)