

Product datasheet for **TS439452P5**

DRP1 (DNM1L) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DNM1L, transcript variant 6, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC239452
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	DRP1 (DNM1L) Mouse Monoclonal Antibody [Clone ID: OTI4F6] (TA506153)
ACCN:	NM_001278465 , NP_001265394
Synonyms:	DLP1; DRP1; DVLP; DYMPLE; EMPF; EMPF1; HDYNIV; OPA5
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001265394
Locus ID:	10059
Cytogenetics:	12p11.21
Protein Pathways:	Endocytosis, Fc gamma R-mediated phagocytosis



[View online »](#)