

## Product datasheet for **TS439766**

### SALL2 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SALL2, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC239766
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_001291446</a></u> , <u><a href="#">NP_001278375</a></u>
Synonyms:	COLB; HSAL2; p150(Sal2); Sal-2; ZNF795
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 $\mu$ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001278375</a></u>
Locus ID:	6297
Cytogenetics:	14q11.2
Protein Families:	Druggable Genome



[View online »](#)