

Product datasheet for **TS439840P5**

NLRX1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NLRX1, transcript variant 5, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC239840
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	NLRX1 Mouse Monoclonal Antibody [Clone ID: OTI4H8] (TA809764)
ACCN:	NM_001282358 , NP_001269287
Synonyms:	CLR11.3; DLNB26; NOD5; NOD9; NOD26
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001269287
Locus ID:	79671
Cytogenetics:	11q23.3
Protein Pathways:	RIG-I-like receptor signaling pathway



[View online »](#)