

## Product datasheet for **TS440090P5**

### AGTPBP1 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of AGTPBP1, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC240090
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	AGTPBP1 Mouse Monoclonal Antibody [Clone ID: OTI9B8] (TA800230)
ACCN:	<a href="#">NM_001286715</a> , <a href="#">NP_001273644</a>
Synonyms:	CCP1; CONDCA; NNA1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001273644</a>
Locus ID:	23287
Cytogenetics:	9q21.33
Protein Families:	Protease



[View online »](#)