

Product datasheet for **TS492563P5**

SARS-CoV-2 N Protein CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of N Protein in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Virus
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone VC102563
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	SARS-CoV-2 N Protein Mouse Monoclonal Antibody [Clone ID: OTI8G10] (TA814432)
ACCN:	<u>YP_009724397</u>
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
Locus ID:	43740575



[View online »](#)