

## Product datasheet for **TS492566**

### S Protein CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of S Protein in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Virus
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone VC102566
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	SARS-CoV-2 S Protein Mouse Monoclonal Antibody [Clone ID: OTI4C6] (TA814463)
ACCN:	<a href="#">YP_009724390</a>
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
Locus ID:	43740568



[View online »](#)