


Product datasheet for TT300006


EZ-LNP Spleen Delivery Kit - 1.0 mL

Product data:

Product Type:	Transfection Reagents
Description:	Instrument-free LNP kit supporting efficient targeted delivery to spleen cells for in vivo applications.
Storage:	4°C for short term (1-2 months); -20°C for long-term. RNA buffer can be stored at 4°C for both short and long-term storage.
Stability:	6 months

Product images:





Add + Mix + Use Workflow.



EZ-LNP Spleen Delivery Kit 1mL.