

Product datasheet for UM570017

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

MEF2C Mouse Monoclonal Antibody [Clone ID: UMAB17]

Product data:

Product Type: Primary Antibodies

Clone Name: UMAB17
Applications: IF, IHC

Recommended Dilution: WB 1:2000, IHC 1:150

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human MEF2C(NP_002388) produced in HEK293T

cell.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.5~1.0 mg/ml (Lot Dependent)

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 51 kDa

Gene Name: myocyte enhancer factor 2C

Database Link: NP 002388

Entrez Gene 17260 MouseEntrez Gene 499497 RatEntrez Gene 4208 Human

006413





Background: This locus encodes a member of the MADS box transcription enhancer factor 2 (MEF2) family

of proteins, which play a role in myogenesis. The encoded protein, MEF2 polypeptide C, has both trans-activating and DNA binding activities. This protein may play a role in maintaining the differentiated state of muscle cells. Mutations and deletions at this locus have been associated with severe mental retardation, stereotypic movements, epilepsy, and cerebral malformation. Alternatively spliced transcript variants have been described. [provided by


RefSeq1

Synonyms: C5DELq14.3; DEL5q14.3


Protein Families: Transcription Factors

Protein Pathways: MAPK signaling pathway

Product images:




Immunohistochemical staining of paraffinembedded Human lymphoma tissue using anti-MEF2C mouse monoclonal antibody. (Clone UMAB17; heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)




Immunohistochemical staining of paraffinembedded Human embryonic brain tissue using anti-MEF2C mouse monoclonal antibody. (Clone UMAB17; heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)






Immunohistochemical staining of paraffinembedded human melanoma using anti-MEFC2 clone UMAB17 mouse monoclonal antibody at 1:200 with Polink2 Broad HRP DAB detection kit ([UM500017]); heat-induced epitope retrieval with citrate buffer, pH6.0 at 95-100C. Nuclear staining seen in the infiltrating lymphocyte and weak staining in tumor cells.




Immunohistochemical staining of paraffinembedded human tonsil using anti-MEFC2 clone UMAB17 mouse monoclonal antibody at 1:200 with Polink2 Broad HRP DAB detection kit ([UM500017]); heat-induced epitope retrieval with citrate buffer, pH6.0 at 95-100C. Nuclear staining seen in both the germinal and non germinal centers of the the tonsil




Immunohistochemical staining of paraffinembedded mouse lung brochus tissue using anti-MEF2C clone UMAB17 mouse monoclonal antibody. HIER TEE buffer pH9 ([B21-100]) at 110C for 10 min, [UM500017] (1:100). Detection was done with Klear Mouse (C/N [D52-18]) DAB Kit.






Immunohistochemical staining of paraffinembedded mouse lung and smooth muscle of blood vessel tissue anti-MEF2C clone UMAB17 mouse monoclonal antibody. HIER TEE buffer pH9 ([B21-100]) at 110C for 10 min, [UM500017] (1:100). Detection was done with Klear Mouse (C/N [D52-18]) DAB Kit.




Immunohistochemical staining of paraffinembedded mouse spleen tissue using anti-MEF2C clone UMAB17 mouse monoclonal antibody. HIER TEE buffer pH9 ([B21-100]) at 110C for 10 min, [UM500017] (1:100). Detection was done with Klear Mouse (C/N [D52-18]) DAB Kit.



Immunofluorescent staining of SH-SY5Y cells using anti-MEF2C mouse monoclonal antibody ([UM500017], green, 1:100). Actin filaments were labeled with Alexa Fluor® 594 Phalloidin (red), and nuclear with DAPI (blue). Scale bar, 20µm.





Confocal immunofluoresce image of primary rat neurons labeled with anti-MEF2C mouse monoclonal antibody ([UM500017], green, 1:50) and with DAPI (blue) for nuclear. Scale bar, 10µm.