

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for UM570080


### ZNF447 (ZSCAN18) Mouse Monoclonal Antibody [Clone ID: UMAB113]

#### **Product data:**


Product Type:	Primary Antibodies
Clone Name:	UMAB113
Applications:	IF, IHC, WB
Recommended Dilution:	IHC 1:200
Reactivity:	Human
Host:	Mouse
lsotype:	lgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human ZSCAN18(NP_076415) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5~1.0 mg/ml (Lot Dependent)
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	54.6 kDa
Gene Name:	zinc finger and SCAN domain containing 18
Database Link:	<u>NP 076415</u>
	Entrez Gene 65982 Human
	<u>Q8TBC5</u>
Synonyms:	ZNF447
Protein Families:	Transcription Factors




## **Product images:**



Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human breast tissue using anti-ZSCAN18 mouse monoclonal antibody. ([UM500080]; heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)




Immunohistochemical staining of paraffinembedded Human Kidney tissue using anti-ZSCAN18 mouse monoclonal antibody. ([UM500080]; heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)




Immunohistochemical staining of paraffinembedded Carcinoma of Human kidney tissue using anti-ZSCAN18 mouse monoclonal antibody. ([UM500080]; heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)

Immunohistochemical staining of paraffinembedded Carcinoma of Human pancreas tissue using anti-ZSCAN18 mouse monoclonal antibody. ([UM500080]; heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)

Immunohistochemical staining of paraffinembedded Human tonsil using anti-ZSCAN18 mouse monoclonal antibody. ([UM500080]; heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)




Western blot of human tissue lysates (15ug) from 10 different tissues (1: Testis, 2: Omentum, 3: Uterus, 4: Breast, 5: Brain, 6: Thyroid, 7: Colon, 8: Spleen 9: Liver, 10: Ovary). Diluation: 1:500.



Immunohistochemical staining of paraffinembedded mouse kidney using mouse anti-ZSCAN18 clone UMAB113 mouse monoclonal antibody. HIER TEE buffer ([B21-100]) (pH9) at 110C for 10 min, [UM500080] (1:100). Detection was done with Klear Mouse (C/N [D52-18]) DAB Chromogen Kit.

Immunohistochemical staining of paraffinembedded mouse small intestine using mouse anti-ZSCAN18 clone UMAB113 mouse monoclonal antibody. HIER TEE buffer ([B21-100]) (pH9) at 110C for 10 min, [UM500080] (1:100). Detection was done with Klear Mouse (C/N [D52-18]) DAB Chromogen Kit.



Immunohistochemical staining of paraffinembedded mouse brain using mouse anti-ZSCAN18 clone UMAB113 mouse monoclonal antibody. HIER TEE buffer ([B21-100]) (pH9) at 110C for 10 min, [UM500080] (1:100). Detection was done with Klear Mouse (C/N [D52-18]) DAB Chromogen Kit.

Immunofluorescent staining of MCF-7 cells using anti-ZSCAN18 mouse monoclonal antibody ([UM500080], green, 1:100). Actin filaments were labeled with Alexa Fluor® 594 Phalloidin (red), and nuclear with DAPI (blue).