


Product datasheet for **UM800179**

HBcAg Mouse Monoclonal Antibody [Clone ID: UMAB289]

Product data:

Product Type:	Primary Antibodies
Clone Name:	UMAB289
Applications:	IHC
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Synthetic peptide (the amino acid sequence is considered to be commercially sensitive) within Human HBcAg(NP_647607). The exact sequence is proprietary.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated

Product images:



Immunohistochemical staining of paraffin-embedded Carcinoma of Human liver tissue using anti-HBcAg mouse monoclonal antibody. (Without high-pressure hot repair, UM800179) (1:3000)



[View online »](#)