

Product datasheet for UM870007

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

MGMT Mouse Monoclonal Antibody [Clone ID: UMAB56]

Product data:

Product Type: Primary Antibodies

Clone Name: UMAB56
Applications: IF, IHC, WB

Recommended Dilution: IHC 1:50~100, IF 1:100

Reactivity: Human, Monkey

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human MGMT (NP_002403) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.5~1.0 mg/ml (Lot Dependent)

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: O-6-methylguanine-DNA methyltransferase

Database Link: NP 002403

Entrez Gene 4255 Human

P16455

Synonyms: methylguanine-DNA methyltransferase; O-6-methylguanine-DNA methyltransferase; O6-


methylguanine-DNA methyltransferase; OTTHUMP0000020741

Protein Families: Druggable Genome






Product images:




Immunohistochemical staining of paraffinembedded Human breast tissue using anti-MGMT mouse monoclonal antibody. ([UM800007]; heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)




Immunohistochemical staining of paraffinembedded Human lung tissue using anti-MGMT mouse monoclonal antibody. ([UM800007]; heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)






Immunohistochemical staining of paraffinembedded Carcinoma of Human lung tissue using anti-MGMT mouse monoclonal antibody. ([UM800007]; heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)




Immunohistochemical staining of paraffinembedded Human thyroid tissue using anti-MGMT mouse monoclonal antibody. ([UM800007]; heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)




Immunohistochemical staining of paraffinembedded Carcinoma of Human thyroid tissue using anti-MGMT mouse monoclonal antibody. ([UM800007]; heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)






Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-MGMT mouse monoclonal antibody. ([UM800007]; heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)




Immunohistochemical staining of paraffinembedded Human lymph node tissue using anti-MGMT mouse monoclonal antibody. ([UM800007]; heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)




Immunohistochemical staining of paraffinembedded Human tonsil using anti-MGMT mouse monoclonal antibody. ([UM800007]; heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)






HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MGMT ([RC229131], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MGMT (1:500).



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-MGMT monoclonal antibody (Clone UMAB56).



Immunofluorescent staining of HeLa cells using anti-MGMT mouse monoclonal antibody ([UM800007], green, 1:50). Actin filaments were labeled with Alexa Fluor® 594 Phalloidin (red), and nuclear with DAPI (blue).