

Product datasheet for UM870021

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Thymidine Phosphorylase (TYMP) Mouse Monoclonal Antibody [Clone ID: UMAB90]

Product data:

Product Type: Primary Antibodies

Clone Name: UMAB90
Applications: IHC, WB
Recommended Dilution: IHC 1:100
Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 311-482 of human

TYMP (NP_001944) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.5~1.0 mg/ml (Lot Dependent)

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 49.8 kDa

Gene Name: thymidine phosphorylase

Database Link: NP 001944

Entrez Gene 1890 Human

<u>P19971</u>

Background: This gene encodes an angiogenic factor which promotes angiogenesis in vivo and stimulates

the in vitro growth of a variety of endothelial cells. It has a highly restricted target cell specificity acting only on endothelial cells. Mutations in this gene have been associated with mitochondrial neurogastrointestinal encephalomyopathy. Multiple alternatively spliced

transcript variants have been identified. [provided by RefSeq, Apr 2012]






Synonyms: ECGF; ECGF1; hPD-ECGF; MEDPS1; MNGIE; MTDPS1; PDECGF; TP

Protein Families: Druggable Genome


Protein Pathways: Bladder cancer, Drug metabolism - other enzymes, Metabolic pathways, Pyrimidine

metabolism

Product images:




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-TYMP mouse monoclonal antibody. ([UM800021]; heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-TYMP mouse monoclonal antibody. ([UM800021]; heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)






Immunohistochemical staining of paraffinembedded Human pancreas tissue using anti-TYMP mouse monoclonal antibody. ([UM800021]; heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)




Immunohistochemical staining of paraffinembedded Human lymphoma tissue using anti-TYMP mouse monoclonal antibody. ([UM800021]; heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)




Immunohistochemical staining of paraffinembedded human tonsil using anti-TYMP clone UMAB90 mouse monoclonal antibody ([UM800021]) at 1:200 with Polink2 Broad HRP DAB detection kit; heat-induced epitope retrieval with GBI Citrate pH6.0 HIER buffer using pre






Immunohistochemical staining of paraffinembedded human breast cancer using anti-TYMP clone UMAB90 mouse monoclonal antibody ([UM800021]) at 1:200 with Polink2 Broad HRP DAB detection kit; heat-induced epitope retrieval with GBI Citrate pH6.0 HIER buffer us





Immunohistochemical staining of paraffinembedded human placenta using anti-TYMP clone UMAB90 mouse monoclonal antibody ([UM800021]) at 1:200 with Polink2 Broad HRP DAB detection kit; heat-induced epitope retrieval with GBI Citrate pH6.0 HIER buffer using p




Immunohistochemical staining of paraffinembedded human lung using anti-TYMP clone UMAB90 mouse monoclonal antibody ([UM800021]) at 1:200 with Polink2 Broad HRP DAB detection kit; heat-induced epitope retrieval with GBI Citrate pH6.0 HIER buffer using press










Western blot of human tissue lysates (15ug) from 10 different tissues (1: Testis, 2: Omentum, 3: Uterus, 4: Breast, 5: Brain, 6: Liver, 7: Ovary, 8: Colon, 9: Spleen, 10: Thyroid). Diluation: 1:500.



Equivalent amounts of cell lysates (10 ug per lane) of wild-type HeLa cells (WT, Cat# LC810HELA) and TYMP-Knockout HeLa cells (KO, Cat# [LC831279]) were separated by SDS-PAGE and immunoblotted with anti-TYMP monoclonal antibody [UM800021] (1:2000`). Then the blotted membrane was stripped and reprobed with anti-HSP90 antibody as a loading control.