

Product datasheet for UM870089

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

Rockville, MD 20850, US
Phone: +1-888-267-4436
https://www.origene.com
techsupport@origene.com
EU: info-de@origene.com
CN: techsupport@origene.cn

PD1 (PDCD1) Mouse Monoclonal Antibody [Clone ID: UMAB197]

Product data:

Product Type: Primary Antibodies

Clone Name: UMAB197
Applications: IF, IHC, WB

Recommended Dilution: IHC 1:200, WB 1:500, IF 1:100...

Reactivity: Human, Mouse, Monkey, Dog, Rat, Rabbit

Host: Mouse Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human PDCD1 (NP_005009) produced in HEK293T

cell

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.5~1.0 mg/ml (Lot Dependent)

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 29.2 kDa

Gene Name: programmed cell death 1

Database Link: NP 005009

Entrez Gene 18566 MouseEntrez Gene 301626 RatEntrez Gene 486213 DogEntrez Gene 5133

<u>Human</u> Q15116





Background: This gene encodes a cell surface membrane protein of the immunoglobulin superfamily. This

protein is expressed in pro-B-cells and is thought to play a role in their differentiation. In mice, expression of this gene is induced in the thymus when anti-CD3 antibodies are injected and large numbers of thymocytes undergo apoptosis. Mice deficient for this gene bred on a BALB/c background developed dilated cardiomyopathy and died from congestive heart failure. These studies suggest that this gene product may also be important in T cell function


and contribute to the prevention of autoimmune diseases. [provided by RefSeq, Jul

Synonyms: CD279; hPD-1; hPD-1; hSLE1; PD-1; PD1; SLEB2


Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Cell adhesion molecules (CAMs), T cell receptor signaling pathway

Product images:




Immunohistochemical staining of paraffinembedded Human tonsil using anti-PDCD1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 3min, [UM800089]) (1:200)




Immunohistochemical staining of paraffinembedded VX2 tumor in rabbit liver using anti-PDCD1 mouse monoclonal antibody. (Heatinduced epitope retrieval by EDTA buffer (pH8.0) at 95-100°C for 30min, [UM800089]) (1:100)






HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PDCD1 ([RC210364], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PDCD1 antibody (1:2000).




Western blot of cell lysates (35ug) from 9 different cell lines (1: HepG2, 2: HeLa, 3: SV-T2, 4: A549. 5: COS7, 6: Jurkat, 7: MDCK, 8: PC-12, 9: MCF7).




Western blot analysis of extracts (35ug) from PBMCs and PHA-stimulated-PBMCs by using anti-PDCD1 monoclonal antibody (1:250).






Western blot of human tissue lysates (15ug) from 10 different tissues (1: Testis, 2: Omentum, 3: Uterus, 4: Breast, 5: Brain, 6: Liver, 7: Ovary, 8: Thyroid, 9: Colon, 10: Spleen). Diluation: 1:500.



Immunofluorescent staining of PDCD1 ([RC210364])-stable-expression cells (left) labeling PDCD1 with mouse monoclonal antibody [UM800089] (1:100, green) and nucleus with Hoechst33342 (blue). HEK293T cells serve as negative control (right).



Immunocytochemistry staining of PDCD1 ([RC210364])-stable-expression cells (left) labeling PDCD1 with mouse anti-PDCD1 monoclonal antibody [UM800089] (1:900). The rihgt is negative control.