

Product datasheet for UM870120

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

PD-L1 (CD274) Mouse Monoclonal Antibody [Clone ID: UMAB228]

Product data:

Product Type: Primary Antibodies

Clone Name: UMAB228
Applications: IHC, WB

Recommended Dilution: WB 1:500, IHC 1:300

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 19-290 of human

CD274 (PD-L1) (NP_054862) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.5~1.0 mg/ml (Lot Dependent)

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 33.28 kDa

Gene Name: CD274 molecule

Database Link: NP 054862

Entrez Gene 29126 Human

Q9NZQ7

Synonyms: B7-H; B7H1; PD-L1; PDCD1L1; PDCD1LG1; PDL1


Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Cell adhesion molecules (CAMs)






Product images:




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human breast tissue tissue using anti-CD274 (PD-L1) mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH9.0) at 120°C for 2.5 min, [UM800120]) (1:300)




Immunohistochemical staining of paraffinembedded Carcinoma of Human liver tissue using anti-CD274 (PD-L1) mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH9.0) at 120°C for 2.5 min, [UM800120]) (1:300)






Immunohistochemical staining of paraffinembedded Carcinoma of Human lung tissue using anti-CD274 (PD-L1) mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH9.0) at 120°C for 2.5 min, [UM800120]) (1:300)




Immunohistochemical staining of paraffinembedded Human lymphoma tissue using anti-CD274 (PD-L1) mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH9.0) at 120°C for 2.5 min, [UM800120]) (1:300)




Immunohistochemical staining of paraffinembedded Human tonsil within the normal limits using anti-CD274 (PD-L1) mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH9.0) at 120°C for 2.5 min, [UM800120]) (1:300)






Immunohistochemical staining of paraffinembedded Carcinoma of Human Esophagus using anti-CD274 (PD-L1) mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH9.0) at 120°C for 2.5 min, [UM800120]) (1:300)




Immunohistochemical staining of paraffinembedded MCF7 (PD-L1 Negative cell line) cell pellets using anti-CD274 (PD-L1) mouse monoclonal antibody (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH9.0) at 120°C for 2.5 min, [UM800120]) (1:300)




Immunohistochemical staining of paraffinembedded NCI-H2228 cell pellets using anti-CD274 (PD-L1) mouse monoclonal antibody (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH9.0) at 120°C for 2.5 min, [UM800120]) (1:300)






Immunohistochemical staining of paraffinembedded MDA-MB-231 cell pellets using anti-CD274 (PD-L1) mouse monoclonal antibody (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH9.0) at 120°C for 2.5 min, [UM800120]) (1:300)




Immunohistochemical staining of paraffinembedded SK-MEL-28 cell pellets using anti-CD274 (PD-L1) mouse monoclonal antibody (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH9.0) at 120°C for 2.5 min, [UM800120]) (1:300)



Immunohistochemical staining of paraffinembedded LoVo (PD-L1 Negative cell line) cell pellets using anti-CD274 (PD-L1) mouse monoclonal antibody (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH9.0) at 120°C for 2.5 min, [UM800120]) (1:300)





Western blot analysis of extracts (35ug) from 5 different cancer cell lines by using anti-CD274 monoclonal antibody (1:500).