

Product datasheet for UM870137

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

AGR2 Mouse Monoclonal Antibody [Clone ID: UMAB245]

Product data:

Product Type: Primary Antibodies

Clone Name: UMAB245

Applications: IHC

Recommended Dilution: IHC 1:500

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 21-175 of human

AGR2(NP 006399) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.5~1.0 mg/ml (Lot Dependent)

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 19.8 kDa

Gene Name: anterior gradient 2, protein disulphide isomerase family member

Database Link: NP 006399

Entrez Gene 23795 MouseEntrez Gene 298961 RatEntrez Gene 10551 Human

095994


Synonyms: AG2; GOB-4; HAG-2; HEL-S-116; PDIA17; XAG-2

Protein Families: Secreted Protein






Product images:




Immunohistochemical staining of paraffinembedded Human breast tissue within the normal limits using anti-AGR2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [UM800137]) (1:500)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-AGR2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [UM800137]) (1:500)






Immunohistochemical staining of paraffinembedded Carcinoma of Human thyroid tissue using anti-AGR2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [UM800137]) (1:500)




Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-AGR2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [UM800137]) (1:500)




Immunohistochemical staining of paraffinembedded Carcinoma of Human bladder tissue using anti-AGR2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [UM800137]) (1:500)






Immunohistochemical staining of paraffinembedded Human Gastric Carcinoma using anti-AGR2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [UM800137]) (1:500).




Immunohistochemical staining of paraffinembedded MCF7 cell pellets using anti-AGR2 mouse monoclonal antibody (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 2.5 min, [UM800137]) (1:500)



Immunohistochemical staining of paraffinembedded T-47D cell pellets using anti-AGR2 mouse monoclonal antibody (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 2.5 min, [UM800137]) (1:500)





Immunohistochemical staining of paraffinembedded MDA-MB-453 cell pellets using anti-AGR2 mouse monoclonal antibody (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.0) at 120°C for 2.5 min, [UM800137]) (1:500)