

Product datasheet for VC100011

pVII (AC_000017) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	pVII (AC_000017) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	pVII
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC100011 represents NCBI reference of AP_000508 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTATATTGATCTCACCTTCTAATAACACTGGATGGGGACTGCGATTCCCCTCCAAGATGTTTGGTG
GCGCAAAAAGAGGAGCGACACCCCGTGAGGGTTCGAGGACATTACCGGGCACCTTGGGGAGCCCA
CAAACGCGGCAGGACAGGACGGACGACTGTGGACGACGCTATCGATGCCGTGGTGAAGAGGCAAGGAAC
TATACTCCCACTCCCCCAGTGTCTACAGTTGACGCTGCCATCCAGACAGTCGTGAGAGGCCACGCA
GCTATGCTAAGATGAAGAGGCGCCGAGGAGGGTGGCCAGGAGGCATAGGAGAAGACCCGGTACCGCGGC
ACAAAGGGCAGCGGCTGCACTGCTCAACCGGGCGCAAGAAGTGGCCGAGGGCCCAATGAGAGCTGCT
AGACGACTGGCAGCCGGCATCGTAACTGTCCCCTCGGAGTAGAAGGAGAGCTGCTGCCGAGCCGCCG
CCGCAATTTCCGCCATGACTCAGGGCAGGCGAGGAAATGTCTACTGGGTGAGAGACTCAGTCTCAGGACT
CAGAGTCCCGGTCAGAACAAGACCACCACGCAAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>VC100011 representing AP_000508 Red=Cloning sites Green=Tags
-------------------	------------------------------------------------------------------

MSILISPSNNTGWGLRFPSKMFGGAKKRSQHPVVRVGRHYRAPWGAHKRGRTGRTTVDDAIDAVVEEARN
YTPTPPPVSTVDAAIQTVVRGARRYAKMKRRRRRVARRHRRRPGTAAQRAAALLNRRARRTGRRAMRAA
RRLAAGIVTVPPRSRRRAAAAAAISAMTQGRGNVYVWRDSVSGLRVVPVTRPPRN

TRTRRLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: AC_000017

ORF Size: 594 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [AC_000017.1](#), [AP_000508](#)

RefSeq ORF: 594 bp

MW: 22.0 kDa