

## Product datasheet for **VC100048**

### 33K (AC\_000007) Virus Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	33K (AC_000007) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	33K
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC100048 represents NCBI reference of AP_000179 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCTCCTAAAAAAGCTGCAGCTCCCCCGCCCCCTCCAACAGATGAGGAGGAGTACTGGGATAGTC  
AGGCTGAAGAGGTGTTGGACGAAGAAGAGGAAATGATGGAGGATTGGGACTCCCTGGACGAGGCCTCTGA  
GGCTGAAGAGGTATCAGACGAAACGCCATCTCCAGCGTGGCCTTCCCATCTCTGCCCCACAGAAACTG  
GCCACGGTACCATCCATAGCTACCACTCCGCCCCACAGGCCCTCCCGCCCTCCCTGTCCAGACGACCAA  
ATCGGCGGTGGGACACGACAGGAAGTACAGCTGCTCCTACTGCACCCGACGAGCTGCTGCAGCTGCTAC  
AGCCGCGTCCACAGAAAGCAGAGAAGACCTGACTCCTAAACCTGACTAAGCCGAAAAAGTCCACCGCG  
GCTGCGGCGGCAGGGGGCGGCGCTGAGACTGGCCCCAAACGAGCCTGTGAGCACTAGAGAGCTGAGAA  
ACCGCATTTTCCAACCTCTACGCAATCTCCAGCAATCTCGGGGTGAGGAACAGGAATTGAAAAATTAA  
GAACCGGAGCCTCGAAGCTTGACCAGATCATGCCTTACCACAAGTCCGAAGACCAGCTGAGGAGAACA  
CTGGAGGACGCTGAGGCTCTTCTCTAAGTACTGCGCCCTGACCCTGAAGGAT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >VC100048 representing AP\_000179  
Red=Cloning sites Green=Tags

MAPKKKLQLPPPPPTDEEEYWDSQAEVLDDEEEMMEDWDSLDEASEAEEVSDETPSPSVAFPPSPAPQKL  
 ATVPSIATTSAPQAPPALPVRPNRRWDTTGTRAAAPTAPAAAAAATAAVTQKQRRPDSKTLTKPKKSTA  
 AAAAGGGALRLAPNEPVSTRELRNRIIPTLYAIFQQSRGQEQELKIKNRSRLRSLTRSCLYHKSEDQLRRT  
 LEDAEALFSKYCALTLKD

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** AC\_000007

**ORF Size:** 684 bp

**OTI Disclaimer:** The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: AC\_000007.1, AP\_000179  
RefSeq ORF: 684 bp  
UniProt ID: P03254  
MW: 25.0 kDa