

Product datasheet for **VC100054**

E1A (AC_000007) Virus Tagged ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | E1A (AC_000007) Virus Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | E1A |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| ORF Nucleotide Sequence: | >The Viral ORF clone VC100054 represents NCBI reference of AP_000161 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s) |

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGACACATTATCTGCCACGGGGTGTCACTGAGGAGATGGCGGCCAGCCTTCTGGATCAGCTGA
TAGAGGAAGTACTTGCCGACAACCTCCCTCCTCACATTTGAGCCACCAACCTTGACGAGCTGTA
CGATCTGGATGTTACGGCCCCGAGGACCCCAACGAAGAAGCCGTATCTCAGATCTTCTGAGTCTGTG
ATGCTCGCTGTGCAAGAAGGCATTGACTTGTTTACGTTCCCTCCAGCCCCGGGTCACCCGAGCCCCCTC
ACTTGAGCCGGCAACCCGAGCAACCAGAGCAGCGCCTTGGGTCCGGTGAGCATGCCAAATCTGGTGCC
AGAGGTATCGATCTGACTTGTACGAGGCCGGCTTCCCCAGCGACGACGAGGATGAGGAAGGCGAA
GAATTTGTGCTGGACTATGTCGAGCACCCCGGGCACGGGTGTCGAGTTGCCATTACCATCGGCGAAACA
CAGGCGACCCAGACATTATGTGTAGCTTGTGTTATATGAGAACTGTGGCATGTTTGTGTATAGCCAGT
GAGTGAACCCGAGCCGAGCCTGAGCCCCGAACCCGAGCCTGCGGACCCACGCGGCGACCAAACTGGTC
CCAGCCATTCTGCGCCGACCAACATCCCCAGTGAAGCCGGAATGTAACAGCAGCACGGATAGTTGTGACT
CAGGTCTAGCAATACCCTCCGAAATCCACCCTGTGGTGCCTTTGCCCTATTAAGCCGTTGCTGT
GCGGGTAGCGGACGAAGACAGGCCGTCGAATGCATTGAAGATCTTCTGAACGAGTCAGGACAACCACTG
GATCTGTCATGCAAAAGACCGAGGCCG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >VC100054 representing AP_000161
Red=Cloning sites Green=Tags

MRHIICHGGVITEEMAASLLDQLIEEVLADNLPPPSHFEPPTLHELVDLVTAPEDPNEEAVSQIFPESV
 MLAVQEGIDLFTFPPAPGSPEPPHLSRQPEQPEQRALGPVSMPLVPEVIDLTCHEAGFPSSDDEDEEGE
 EFVLDYVEHPGHGCRSCHYHRRNTGDPDIMCSLCYMRTCGMFVYSPVSEPEPEPEPEPARPTRRPKLV
 PAILRRPTSPVSRECNSTSDSCSGPSNTTPEIHPVVPLCPKIPVAVRVGGRRQAVECIEDLLNESGQPL
 DLSCKRPRP

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: AC_000007

ORF Size: 867 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: AC_000007.1, AP_000161
RefSeq ORF: 867 bp
UniProt ID: P03254
MW: 31.9 kDa