

Product datasheet for **VC100055**

E1B 19K (AC_000007) Virus Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: E1B 19K (AC_000007) Virus Tagged ORF Clone
Tag: Myc-DDK
Symbol: E1B 19K
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >The Viral ORF clone VC100055 represents NCBI reference of AP_000162 with codon optimized for human cell expression
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**

ATGGAGGCCTGGGAATGTCTGGAGGACTTCAGCGCTGTCCGGAATCTGCTTGAACAGAGCAGCAATAGTACTTCATGGTTCTGGAGGTTTCTGTGGGGTCCAGTCAAGCGAAGCTTGTCTGCCGAATCAAAGAAGACTATAATGGGAATTTGAAGAGCTGCTCAAGTCATCGGGGAGCTGTTTCGATAGTCTGAATCTCGGACACCAGGCTCTCTCCAAGAGAAGGTGATTAACACTCTGGATTTTTCAACACCAGGACGCGCAGCGCCCGCTGGCTTTTCTGAGCTTCATAAAAGATAAGTGGAGCGAAGAACTCACTTGTGTCAGGTGGCTATCTGCTGGATTTCTGGCGATGCATCTGTGGCGGAGTGGTACGCCACAAGAACCAGCTCCTCCTGCTGTCTAGTGTGCGGCCCGCCATCATCCCTACCGAGGAACAGCAGCAGGAAGAGGCCAGGAGAAGACGACGACAAGAACAGAGTCCCTGGAACCCCGGGCAGGTCTTGACCCTAGAGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >VC100055 representing AP_000162
Red=Cloning sites **Green**=Tags

MEAWECLEDFSAVRNLLQSSNSTSWFWRFLWGSSQAKLVCRKEDYKWEFEELLKSCGELFDSLNLGHQALFQEKVIKTLDFSTPGRAAAVAFLSF IKDKWSEETHLSGGYLLDFLAMHLWRAVVRHKNRLLLLSSVRPAIIPTEEQQQEARRRRRQEQSPWNPRAGLDPRE

TRTRRLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:

Cloning sites used for ORF Shutting:



* The last codon before the Stop codon of the ORF

ACCN: AC_000007

ORF Size: 525 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [AC_000007.1](#), [AP_000162](#)
RefSeq ORF: 525 bp

UniProt ID: [P03254](#)
MW: 20.5 kDa